

OS基準のLabchart, Chart対応表

		Labchart v7	Labchart v6	Chart v5	Chart v4	Chart v3
Windows	windows 7	v7.0.3~				
	windows Vista	○	○ *1	v5.5.1~		
	windows XP	○ (SP3~)	○ (SP2~)	○	v4.1.2~	
	windows 2000			○	○	○
	windows Me				○	○
	windows NT4				○	○
	windows 98				○	○
	windows 95				○	○
Mac	Mac OSX 10.7~	v7.2.3~				
	Mac OSX 10.5, 10.6	○	○	v5.5.6		
	Mac OSX 10.4~	○	○	○		
	Mac OSX 10.1~			○		
	Mac OSX Classic環境				○	
	Mac OS8.6~9.2				○ (PowerPC以降)	○ (86k以降)

*1 v6.1.1~ 64-bitにも対応